




venome®



VENOME PEN

MICRONEEDLING PEN

MICRO-NEEDLE MESOTHERAPY

One of the most effective methods to achieve strong skin stimulation and revitalisation, involving controlled microperforation of the epidermis without the risk of damaging blood vessels.

The biggest advantage of micro-needle mesotherapy is the ability to control the intensity and depth of the treatment in the same way, across the entire treatment area. This type of treatment is an ideal complement to needle mesotherapy treatments, enhancing the final result.



MICRO-NEEDLE MESOTHERAPY MECHANISM

Controlled micro-damage to the skin through targeted puncturing

Release of blood droplets from blood vessels initiates a sequence of thrombocyte activation processes

Thrombocytes release growth factors

Growth factors intensify and stimulate restoration and complete remodelling processes by stimulating type III collagen production

INDICATIONS

Wrinkles

Uneven skin texture

Flaccidity

Discolouration

Excessive dryness

Stretch marks

Scars

NEEDLE LENGTH

0.20 - 0.49 mm

MECHANISM

Increased penetration of active substances

USE

Enhanced results of products used

0.50 - 0.90 mm

Increased absorption of active substances and stimulation of collagen production

Pigmentation, fine lines, hair loss

1.00 - 2.00 mm

Increased penetration of active substances and natural collagen induction

Wrinkles, scars, acne, pigmentation, lifting, smoothing out texture

VENOME PEN

Venome Pen Ultima A6-S is a device designed to perform micro-needle mesotherapy treatments with the use of carefully selected cartridges and adjustable skin perforation depth.

Needle extension range: 0.25 - 2.5 mm

Puncture speed: 22,000 punctures per minute

Wireless working time up to 8 hours

Venome Pen is the perfect tool in the fight against wrinkles, acne scars, pigmentation or stretch marks. The device can operate in both wired and wireless mode, offering great flexibility.

HOW DOES THE **VENOME PEN** WORK?

Micro-puncture stimulates the skin's natural collagen and elastin production processes, making the skin firm and supple. Thanks to the punctures made, nutrients from applied product reach deeper into the skin, effectively nourishing it.

ADVANTAGES OF USING MICRO-NEEDLE MESOTHERAPY WITH THE VENOME PEN

Quick treatment

Precise treatment

Injection depth control

Possible to make more punctures

Intense stimulation of the skin on the same levels

ADVANTAGES OF THE VENOME PEN

High puncture rate of up to **22,000 punctures** per minute

LCD display indicating the current operating mode (1-6)

High-quality aluminium casing

Wired and wireless operation mode

Stable 3-point cartridge attachment

Cartridges with the thinnest needles of 0.18 mm - making the treatment more effective and less painful

SET INCLUDES

- Venome Pen Ultima A6-S device
- 2 pcs of Lithium-ion batteries (500 mah)
- Wired operation module
- Charger
- 2 pcs of 16PIN cartridges
- Power supply + cable

Product intended for professional use.