




HAIR? WHERE? *EVERYWHERE!*

INTENSIVE HAIR LOSS TREATMENT

INDICATIONS



WEAKENED DISULPHIDE BRIDGES IN THE HAIR



WEAK AND FRAGILE SHAFTS PRONE TO BREAKING



STRESS- OR MEDICATION-RELATED HAIR LOSS



EXCESSIVE LOSS



WEAK ROOTS AND SHAFTS AFTER HAIR TREATMENTS



ANDROGENETIC PROLAPSE

PRODUCTS



VENOME MESO HAIR



VENOME HAIR FILLER 2.0


A **trichological therapy** with revitalising properties, combining Venome Hair Filler 2.0 and Venome Meso Hair as a unique treatment for weak hair prone to excessive loss.

Venome Hair Filler 2.0 is a unique combination of hyaluronic acid, hydrolysed keratin and growth factors, and biomimetic peptides that effectively addresses the problem of falling out and weakened hair, stimulating the revitalisation of hair roots and the regrowth of healthy hair. The growth factors and peptides used in the formula work synergistically at different levels, stimulating cells to grow and work continuously, respectively.

Venome Meso Hair is the ideal formulation for strengthening weak and fragile hair. Thanks to its highly concentrated formula of keratin, biotin and a composition of vitamins B6, B1 and panthenol, the mesococktail supplements vitamin deficiencies, improves hydration, strengthens hair elasticity and minimises hair loss.

COURSE OF TREATMENT

1. **USING A COLLECTION NEEDLE**, TAKE 2ML OF VENOME HAIR FILLER 2.0 INTO A 5ML SYRINGE.
2. **USING A COLLECTION NEEDLE**, TAKE 2ML OF VENOME MESO HAIR INTO A NEW 5ML SYRINGE.
3. **CONNECT** THE TWO SYRINGES WITH **A CONNECTOR**.
4. **WITH A GENTLE MOVEMENT** OF THE SYRINGE PLUNGER, **MIX THE TWO PREPARATIONS TOGETHER**.
5. **TRANSFER THE FINISHED PRODUCT** INTO ONE OF THE SYRINGES.
6. **REMOVE THE CONNECTOR**.
7. **DISINFECT** THE TREATMENT AREA.
8. **PROCEED WITH A MESOTHERAPY TREATMENT**.



POST-TREATMENT CARE

AFTER THE TREATMENT, **DO NOT TOUCH THE TREATED AREAS** FOR FASTER RECOVERY.